

Excellence in
Infection Control



Bacillocid[®] Extra

Modern aldehyde containing surface disinfectant. **Bacillocid[®] Extra** has an excellent cleansing property, and a broad spectrum of germicidal action.



Effective against HBV, HIV, HCV, AVIAN FLU, SWINE FLU, SARS, SPORES
Passes EN 13727, EN 13624, EN 14348, EN 1040, EN 1275, EN 14347, EN 14476



Bacillocid® Extra

Modern aldehyde containing surface disinfectant. **Bacillocid® Extra** has an excellent cleansing property, and a broad spectrum of germicidal action.

Recommended by:



Product Properties

short exposure times

Effective against

Clostridium Difficile Spores

Excellent cleaning power

Composition

Each 100gm contains :

(Ethylenedioxy) dimethanol - 14.1g

Glutaraldehyde - 5.0g

Specialized, Hi-Tech, Inbuilt

Cleansors

Corrosion Inhibitors

Microbiology

Bactericidal, Yeasticidal

Fungicidal, Tuberculocidal

Mycobactericidal, Sporicidal

Virucidal against enveloped

viruses (incl. HBV, HIV, HCV),

Virucidal

Areas of application

Bacillocid® Extra is suitable for the cleaning surface disinfection of washable surfaces using the wet-wipe procedure, for example : for medical devices and inventory, which come under the Medical Device Directive in hospitals and residential / nursing homes, particularly for all functional areas (acc. to BPD*) for areas in the pharmaceutical industry that are relevant to hygiene (acc. to BPD*) in laboratories and the cosmetics industry (acc. to BPD*)

Directions for use

Bacillocid® Extra is supplied as a concentrate. Prepare use-solutions with cold water. Completely wet the washable surfaces (e.g. floors) with an adequate amount of solution.

Material compatibility

Bacillocid® Extra use-solutions were tested for their compatibility on the following materials among others :

Metals : Stainless steel (V2A), aluminium, copper, brass.

Plastics : Polyethylene, polypropylene, polystyrene, polyurethane, PVC, silicon, rubber latex, Makrolon®, acrylic glass, Teflon®.

Proven efficacy

Bacteria and Fungi	
EN Phase 2 / Step 1 Efficacy according to EN Norm Phase 2 / Step 1 (suspension tests), tested under clean / dirty conditions	Bactericidal (EN 13727) - clean conditions - dirty conditions Yeasticidal (EN 13624) - clean conditions - dirty conditions Fungicidal (EN 13624) - clean conditions - dirty conditions Tuberculocidal (EN 14348) - clean conditions Mycobactericidal (EN 14348) - clean conditions
EN Phase 1 / Basic tests Efficacy according to EN Phase 1 (Basic tests / suspension tests) without contamination; does not define the applicability of a product for a specific purpose	Bactericidal (EN 1040) Yeasticidal (EN 1275) Fungicidal (EN 1275)
VAH Certified application recommendations for prophylactic wet-wipe disinfection from the Association for Applied Hygiene (VAH). Based on suspension and practical tests, tested under clean conditions (i.e. optically clean surfaces)/ dirty conditions (i.e. visibly contaminated surfaces)	Bactericidal / Yeasticidal - clean conditions - dirty conditions Fungicidal - clean conditions - dirty conditions Tuberculocidal - clean conditions - dirty conditions
RKI Recognized substance for decontamination acc. to § 18 IfSG (Robert Koch-Institute [RKI]) Area B - see 'Viruses'	Area A - veg. Bacteria (incl. Mycobacteria), fungi, and fungi spores
Bacterial spores	
Efficacy against spores	Sporicidal (EN 14347) C.diff-spores (Ribotype 027)
Viruses	
EN Phase 2 / Step 1 Efficacy according to EN Phase 2 / Step 1 (suspension tests), tested under clean / dirty conditions	Virucidal (EN 14476) - clean conditions - dirty conditions Adeno virus (EN 14476) - clean conditions - dirty conditions
Efficacy against viruses (German Society for the Control of Viral Diseases (DVV))	Virucidal against enveloped viruses (incl. HBV, HIV, HCV) Virucidal
Appraised efficacy against non- enveloped viruses (in accordance with (DVV))	Adeno virus Polyoma virus
Appraised efficacy against non- enveloped viruses (in accordance with (DVV))	HAV Rota virus
Appraised efficacy against non- enveloped viruses (in accordance with EN)	MNV - clean conditions - dirty conditions
Recognized substance for decontamination acc. to § 18 IfSG (RKI) Area A - see 'Bacteria'	Area B - enveloped and non - enveloped viruses

Dilution :

**Wet Wipe, Mopping / Terminal
Disinfection Non Critical Areas**

0.12.5% (Clean) = 1.25ml x 1 litre

0.25% (Dirty) = 2.5ml x 1 litre

**Wet Wipe, Mopping / Terminal
Disinfection Critical Areas**

0.50% (Clean) = 5ml x 1 litre

1% -2% (Dirty) = 10ml - 20ml x 1 litre

Fogging of Critical Areas

0.50% (Clean) = 5ml x 1 litre

1% -2% (Dirty) = 10ml - 20ml x 1 litre

**Mopping (Wet Wipe) is
recommended over fogging.**

CDC - USA - Guidelines: 2003, 2008

"Do not perform disinfectant fogging for
routine purposes in patient-care areas."
Category IB, Category II.



Packing :

Content

500ml

5 Ltrs

Carton

20 bottles

2 Jars



Application times for final disinfection

After discharge / transfer of patients infected and colonised by pathogens respectively, the patient room needs to be disinfected. The final disinfection is carried out in areas or rooms that have been used for patient care or treatment and has to include all surfaces and objects near the patient that are or could be contaminated (1). In case of epidemics, the selected product needs to be listed by the Robert Koch-Institute (RKI) - for norovirus outbreaks, for instance, the final disinfection has to be performed with a virucidal disinfectant listed for activity areas A (bacteria) and B (viruses) (2).

1. Empfehlung der Kommission für Krankenhaushygiene und Infektionsprävention beim Robert Koch-Institut (RKI). Anforderungen an die Hygiene bei der Reinigung und Desinfektion von Flächen. Bundesgesundheitsbl - Gesundheitsforsch - Gesundheitsschutz 2004 47:51-61
2. Liste der vom Robert Koch-Institut gerprüften und anerkannten Desinfektionsmittel und -verfahren. Bundesgesundheitsbl - Gesundheitsforsch - Gesundheitsschutz 2004 : 47:51-61.2 Stand vom 31.5.2007 (15. Ausgabe)

Exclusively Manufactured & Marketed by:

RAMAN & WEIL PVT. LTD.

Registered Office : 15, Chateau Marine, B-Rd., Churchgate, Mumbai - 20.

Administrative Office : Room No. 36/37, 3rd Floor, Mittal Chamber,

Opp. Inox Theatre, Nariman Point, Mumbai - 400 021, India

Phone : 22049527, 22856397 Email : sales@ramanweil.com

Visit : www.ramanweil.com



In technical collaboration with:

BODE Chemie GmbH

Melanchthonstraße 27

D-22525 Hamburg, Germany

Must Visit : www.bode-science-center.com

Recipient : ISO 9001, ISO 14001, ISO 13485

EC Certificate of Conformity, WHO - GMP